




MR Solutions for MR RT

Klarity MR Solutions consists of MR-Safe Overlays and Coil Supports which are designed to enable MR imaging in the same position as treatment.




FEATURES:

- Klarity's high quality MR Overlays provides strong support and an even surface.
- The Overlay is easy to install and accommodates a variety of MR-compatible positioning and immobilization devices.
- The MR-Safe Coil Support is height-adjustable and can be adjusted to bring the coil close to each individual patient.

Customization is available from Klarity for your specific needs.

KLARITY MR SOLUTION FOR PHILIPS INGENIA 3.0T MR



Solution Code: R625-MPGF

Components & Specifications:

MR OVERLAY (1)

- Dimensions: 213x53x2.2cm
- Weight: 7.0 kg
- Material: Fiberglass + Kevlar
- The Overlay accommodates a variety of positioning and immobilization devices.

BODY COIL SUPPORT (1)

- Material: ABS + Fiberglass
- Height adjustable and can be adjusted to bring the coil close to each individual patient.
- The Support is easy to fix on the overlay and holds the coil inside without compromising the image quality.

MR INDEXING BAR (2)

- Item Number: RE-7BGF

KLARITY MR SOLUTION FOR SIEMENS MAGNETOM SKYRA



Solution Code: R625-MSGF

Components & Specifications:

MR OVERLAY (1)

- Dimensions: 205.5x58x2.2cm
- Weight: 6.6 kg
- Material: Fiberglass + Kevlar
- The Overlay accommodates a variety of positioning and immobilization devices.

HEAD COIL SUPPORT (1)

BODY COIL SUPPORT (1)

- Material: ABS + Fiberglass
- Height adjustable and can be adjusted to bring the coil close to each individual patient
- The Support is easy to fix on the overlay and holds the coil inside without compromising the image quality.

MR INDEXING BAR (2)

- Item Number: RE-630EGFS