

Sheets and Bolus

Klarity® Thermoplastic Sheets

Klarity low-temperature thermoplastic material is designed specifically for radiation therapy. It softens in warm water or in a dry-heat oven and molds quickly to facilitate treatment. It provides sufficient setup time for the therapist, and simply reheating the thermoplastic will allow for any necessary modifications.

Klarity® EZ Bolus™ Thermoplastic Bolus in sheets and pellets

Available in sheet and pellet form, Klarity EZ Bolus™ is designed as an aid to spare skin and treat irregular contours or locations. After heating to about 150°F (65°C), EZ Bolus™ is easy to manipulate and shape by hand, making it ideal for eliminating air gaps. Dosimetric properties are the same as conventional mask thermoplastics.

R-4000 Klarity EZ Bolus, solid, 4.0 mm, 92 x 61 cm

R-4000A Klarity EZ Bolus, solid, 4.0 mm, 61 x 46 cm

R-4000E Klarity EZ Bolus, solid, 4.0 mm, 15 x 15 cm


R-4000F Klarity EZ Bolus, solid, 4.0 mm, 30 x 30 cm

R-4000G Klarity EZ Bolus, solid, 4.0 mm, 20 x 23 cm

R-601 Klarity EZ Bolus Thermoplastic Pellets, 1.1-lb jar

R-603 Klarity EZ Bolus Thermoplastic Pellets, 3.3-lb jar

PERFORATIONS:



15% perf

36% perf

22% perf

42% perf


THICKNESSES:

1.6 mm
(1/16 in)

2.0 mm
(1/12 in)

2.4 mm
(3/32 in)

3.2 mm
(1/8 in)



Sheets and Bolus

Bolx I™ and Bolx II™

Bolx I and Bolx II are flexible build-up tissue equivalents. Both are designed to increase the surface dose, and both offer excellent conformity to body contours.


Both Bolx I and II are made of a cured transparent gel. Bolx I is encased in a layer of polyethylene skin. Bolx II is skinless and is intended for single-patient use.

Bolx I™ Bolus Material (with polyethylene skin)

- KAC-31052** Bolx I Sheet, 30 x 30 x 0.5 cm
- KAC-31102** Bolx I Sheet, 30 x 30 x 1 cm
- KAC-31152** Bolx I Sheet, 30 x 30 x 1.5 cm
- KAC-31202** Bolx I Sheet, 30 x 30 x 2 cm
- KAC-31302** Bolx I Sheet, 30 x 30 x 3 cm

Bolx II™ Bolus Material (without skin)

- KAC-32052** Bolx II Sheet, 30 x 30 x 0.5 cm
- KAC-32102** Bolx II Sheet, 30 x 30 x 1 cm
- KAC-32152** Bolx II Sheet, 30 x 30 x 1.5 cm
- KAC-32202** Bolx II Sheet, 30 x 30 x 2 cm
- KAC-32302** Bolx II Sheet, 30 x 30 x 3 cm



Gel Bolus with removable skin

This cost-effective bolus is made of transparent gel which comes with removable polyethylene skin. So it can be used with or without its skin.

The bolus is flexible build-up tissue equivalent designed to increase the surface dose and offer excellent conformity to body contours.

The gel bolus is intended for single-patient use.



- RBL-T2815** 28 x 28 x 1.5 cm
- RBL-T2810** 28 x 28 x 1.0 cm
- RBL-T2808** 28 x 28 x 0.8 cm
- RBL-T2805** 28 x 28 x 0.5 cm
- RBL-T2803** 28 x 28 x 0.3 cm
- RBL-T2015** 20 x 20 x 1.5 cm
- RBL-T2010** 20 x 20 x 1.0 cm
- RBL-T2008** 20 x 20 x 0.8 cm
- RBL-T2005** 20 x 20 x 0.5 cm
- RBL-T2003** 20 x 20 x 0.3 cm