

Klarity[®] Optek™ System


Optimal stability for the head, neck, and shoulders

- Enhanced indexing
- Comfortable shoulder retraction
- Strong and durable



Klarity Optek System R630-3SCF

The Klarity Optek System is a comfortable and adaptable S-type Overlay for head and neck treatments. The Klarity-K™ Carbon Fiber construction has enhanced strength and durability. The system is ideal for IMRT and conformal treatments in a 360° arc and is dual-configured for adult and pediatric treatments. The lower portion gently tapers for patient comfort and the overlay is equipped with additional indexing around the edges. The Optek™ System is extremely lightweight and CT-compatible with excellent attenuation factors. Comfortable and adjustable hand grips assist with shoulder retraction. The Optek™ System is compatible with all S-type masks.

105cm (Length) x 61cm (Width)
Also available as MR Compatible 

INCLUDED COMPONENTS:



Overlay



Plastic Handgrips R630-TPM

OPTIONAL COMPONENTS:



**Carbon Fiber Handgrips
R630-TPCF**



**Aluminum Indexing Bar
RE-630EA**



**Prone Head MultiFix™
R605-HPCF**



**Tilting Head MultiFix™
R602-DCF**

COMPATIBLE WITH:



S-type Head Masks



S-type H&S Masks



**Indexing Bars
RE-630EA , RE-7E**



**Head Supports
R507 Series**





Klarity AccuCushions®



Vacuum Bags

Klarity[®] Optek Demi-Baseplate

The Klarity Optek™ Demi-Baseplate for the Elekta HexaPOD™ System.



Klarity Optek Demi-Baseplate

The Optek Demi-Baseplate is specially designed for the Elekta HexaPOD™ evo RT System, and is equally useful for situations requiring a shortened baseplate. It is extremely lightweight and the lower portion gently tapers for patient comfort.

The Optek™ Demi-Baseplate is available with 3 types in order to accommodate different types of masks.

R630-SHCF, compatible with S-Type masks

R630-VHCF, compatible with V-Type masks

R630-OHCF, compatible with O-Type masks

58cm (Length) x 61cm (Width)



COMPATIBLE WITH:



Indexing Bars



Head Supports



Klarity AccuCushions[®]



Vacuum Bags