

Klarity®

Klarity MR Vacuum Bags

- Fully MR-safe
- Airtight and ultra-rigid
- A wide selection of shapes and dimensions, customized service is offered



Ordering Information of Klarity MR Vacuum Bags

Model	Size	Suitable
MR7202-9NM	Vacuum bag for Head & Shoulder Support, 65.2X68.5CM	Universal
MR7206-16NM	Vacuum bag for Head & Shoulder Support, 98.5X70CM	Universal
MR7244-21NMR-O	Vacuum bag for Head & Shoulder Support, 135.2X60CM	for R612-NCF
MR7507-36NM	Vacuum bag for Body support, 123X81.5CM	Universal
MR7504-28NM	Vacuum bag for Body support, 100X70CM	Universal
MR7580-36NM	Vacuum bag for Body support, 140X70CM	Universal
MR75103-77NM	Vacuum bag for Body support, 200X80CM	Universal
MR7581-80NMA	Vacuum bag for Body support, 200X88CM	Compatible with Elekta bodyfix
MR7512-45NM	Vacuum bag for Body support, 150X80CM	Universal
MR7644-40NM	WingSpan™ Vacuum Bag, 142X92CM	Universal
MR7641-25NMB2	WingSpan™ Vacuum Bag, 120X80CM	Compatible with Klarity SBRT System
MR7704-5NMT-O	Vacuum bag for Body support, 80X28CM	Compatible with Klarity SBRT System

* Available with valves to fit all standard vacuum pumps



Valves(R7800-8)



Durable nylon shell material



Low-density mini-spheres polystyrene



Non-magnetic valves

