

Klarity®

Klarity AXIS Patient Setup Platform





Klarity AXIS Patient Setup Platform KT-811B



- **Improves Work Efficiency and Patient Throughput**

Klarity AXIS is an independent patient setup platform which improves radiotherapy centers' efficiencies and their patient throughput



- **Human-centered Design**

1. Sensitive touch screen allows operation when wearing medical gloves
2. Electric lifting height of the treatment couch ranges from 65cm to 100cm
3. Couch Dimensions: 220cm (Length) x 53cm (Width)



- **The Treatment Couch is Translational in Three Axes**

1. Moving forwards and backwards along the X axis (surging) from 0 to 70cm
2. Moving left and right along the Z axis (swaying) from -20 to 20cm
3. Moving up and down along the Y axis (heaving) from 65 to 100cm

Klarity Plain Patient Setup Platform KT-811A

- Couch Dimensions: 200 x 53 x 80cm
- Come with a step ladder for patient's easy getting on and off the couch

