

Klarity

WingSpan™ Baseplate



- Secure arm support
- Lightweight and easily transportable
- Low attenuation

Klarity WingSpan™ Baseplate R632-BTCF

The Klarity WingSpan™ Baseplate features secure arm support for upper body treatment and diagnostic procedures. The WingSpan™ Baseplate comes with T-Hand Grips. It is compatible with standard Silverman head supports and can be used with the Head Cup included with the Supine Breast System. The compatible WingSpan™ Vacuum Bag allows for secure customizable support. It is comfortable and lightweight, and allows for easy transporting, indexing, and re-positioning, with low attenuation. The Klarity WingSpan™ Baseplate can be used alone, or in conjunction with the Klarity Supine Breast System. It is also compatible with the Klarity SBRT System.

56.1cm (Length) x 58cm (Width)

Also available as MR Compatible 

COMPATIBLE WITH:



Supine Breast System



SBRT System



Head Supports



Indexing Bars



Vacuum Bags