

Klarity[®]

Carbon Fiber CT Overlays

- Made of high quality materials with a high-pressure vacuum forming process
- Uniform density inside avoids unacceptable artifacts during CT scanning
- Easy to install with quick-lock mechanism
- Stable, reliable, and flexible to adjust for levelling



Customized service for CT Overlays is available.

CT Brand	CT Model	Klarity Overlay	Spec. (mm)
Siemens	Somatom Emotion	R625-SFCF	2140x530x35
	Somatom Scope		
	Somatom Sensation	R625-EFCF	2160x530x40
	Somatom Definition AS	R625-SFCF3	2140x530x35
	Somatom Definition Flash	R625-SFCF4	2140x530x35
	Somatom go.Top	R625-SFCF5	2140x530x35
Philips	MX 16 - slice	R625-PFCF	2140x530x35
	Brilliance Big Bore	R625-PFCF1	2500x530x25
	Brilliance 6 Slice CT	R625-PFCF2	2500x530x25
Neusoft	NeuViz 64 IN	R625-NFCF	2140x530x35
Toshiba	Aquilion one	R625-NFCF2	2140x530x35
United Imaging	uCT 550	R625-PFCF3-A1	2140x530x35
GE	GoldSeal Optima CT660	R625-GFCF1	2140x530x35
Hitachi	Supria	R625-HFCF	2140x530x35

