

Klarity[®]

PROMISE[™]

Klarity PROMISE[™] Proton Therapy System





Klarity PROMISE™ Proton Therapy System



PROMISE™ Baseplate
(Available with versions
for CT sim, MR sim and
treatment delivery)



PROMISE™ Head AccuCushion
Item No.: RG-550-CEL2



PROMISE™ Head Mask
Item No.: RG-PH-2436C

Klarity PROMISE™ Proton Therapy Solutions

■ For CT sim



Components:

1. PROMISE™ Baseplate for CT sim: R680-ACF
2. PROMISE™ Head Mask: RG-PH-2436C
PROMISE™ H&S Mask: RG-PS-2436C
3. PROMISE™ Head AccuCushion: RG550-CEL2
4. PROMISE™ Indexing Bars: RE-7S, RE-630ES

Features:

- Material: Carbon Fiber
- Dimensions: 1203x530x22mm
- The mask can be easily and steadily fixed with the baseplate
- Adjustable hand-grips enable easier and personalized immobilization
- Can be fixed onto the couchtop by two indexing bars for best stability

■ For MR sim



Components:

1. PROMISE™ Baseplate for MR sim: R680-AGF
2. PROMISE™ Head Mask: RG-PH-2436C
PROMISE™ H&S Mask: RG-PS-2436C
3. PROMISE™ Head AccuCushion: RG550-CEL2
4. PROMISE™ Indexing Bars: RE-7SGF, RE-630ESGF

Features:

- Material: Fiber Glass
- Dimensions: 1203x530x22mm
- The mask can be easily and steadily fixed with the baseplate
- Adjustable hand-grips enable easier and personalized immobilization
- Can be fixed onto the couchtop by two indexing bars for best stability

■ For treatment delivery



Components:

1. PROMISE™ Baseplate for proton treatment delivery: R680-ACFP
2. PROMISE™ Head Mask: RG-PH-2436C
PROMISE™ H&S Mask: RG-PS-2436C
3. PROMISE™ Head AccuCushion: RG550-CEL2

Features:

- Material: Carbon Fiber
- Dimensions: 1577x530x28mm
- The mask can be easily and steadily fixed with the baseplate
- Adjustable hand-grips enable easier and personalized immobilization
- Baseplate connected with and locked into the couchtop with two metallic pins



Klarity Medical & Equipment (GZ) Co., Ltd

No.14, 3rd Street Shawan, GETDD Guangzhou 510730, China

Tel: +86-20-8222 6380

Fax: +86-20-8208 8552

Email: sales@klarity-medical.com

www.Klarity-Medical.com