




# MR Patient Positioning Devices





# MR Patient Positioning Devices

## Baseplates





## MR Vacuum Bags

The newly developed Klarity MR Vacuum Bags are made of MR safe fabrics and filled with MR safe micro-balls of Styrofoam. To ensure its safe use, in the same time to guarantee the imaging quality, we equip each bag with a special valve which is MR conditional up to 3.0T.



## Accessories





**MR Movable KneeFix Set**  
R612-LGF

**MR AIO Baseplate**  
R612-NGF

**MR Overhead Arm Positioner**  
R612-WGF




## Klarity MR Solo Align™ Full Body System

- An all-in-one baseplate for S-type masks compatible with small-bore machines.
- A/B/C location choices to provide better immobilization for patients.
- The MR compatible baseplate can be used to immobilize the head, neck, shoulders, thorax, and pelvis for targeted as well as whole-body treatments.

**MR**  
**R612**

### System Features



An adjustable head frame to fit different body sizes.



A groin lock enhances pelvic mask fixation, which can be easily mounted/dismounted.



Three pairs of indexing holes.





# One product, many solutions



### For Head & Neck Fixation

**Required**

Baseplate R612-NGF  
Head Mask: RLG460SHT (LiteGreen™)  
Narrow Shoulder H&S Mask:  
RLG460-2436S5T (LiteGreen™)



### For Head & Thorax Fixation

**Required**

Baseplate R612-NGF  
Overhead Arm Positioner R612-WGF  
Head & Thorax Mask :  
RLG322-S42436 (LiteGreen™)

\*A supine positioning solution for radical mastectomy IMRT




### For Breast Fixation

**Required**

Baseplate R612-NGF  
Overhead Arm Positioner R612-WGF  
Mask R322-S424C-L (Left Open)  
or R322-S424C-R (Right Open)

\*A supine positioning solution for breast-conserving IMRT



### For Craniospinal Irradiation

**Required**

Baseplate R612-NGF  
Narrow Shoulder H&S Mask:  
RLG460-2436S5T (LiteGreen™)  
Pelvis Mask RLG902-2436BA (LiteGreen™)



### For Pelvis Fixation

**Required**

Baseplate R612-NGF  
Overhead Arm Positioner R612-WGF  
Pelvis Mask: RLG902-2436BA (LiteGreen™)



### For Multi-Target RT Body Fixation

**Required**

Baseplate R612-NGF  
Overhead Arm Positioner R612-WGF  
Trunk Mask: RLG902-2436LA (LiteGreen™)

MR



### Klarity MR SRS System

- Double-masks system molds to the front and back of the head, enhancing accuracy and reliability.
- Frameless, non-invasive patient positioning system with sub-millimeter accuracy.
- The unique SRS Shim System is very easy to use with three alternative frames

MR  
SRS

#### • Components



Klarity MR SRS Baseplate  
R660-AGF1



Alternative Frames, a Set of 3 types.

#### SRS Consumables Package



SRS Masks


Code: RBSRS-02

+



BiteLok™

+



SRS Openface Masks

Code: RBSRS-F02



MR Patient  
Positioning Devices



MR  
SBRT

### Klarity MR SBRT System

The Klarity MR SBRT System is simple to use, highly indexable, and easily adjustable. The baseplate, knee and leg positioner, and vacuum bags are dedicated to maximize patient comfort during treatment.

- **Components**



Belly Compression Paddle and  
Respiratory Belt for Breath Restriction

WingSpan™ Baseplate

Foot Holder

Vacuum Bag for Body Support

Vacuum Bag for Knee Support

**Klarity MR SBRT System**

Product Code: R62400-GF

(Please Specify When You Choose Klarity  
SBRT System for Philips MR Sim)



**Klarity PROMISE™ Proton Therapy System**

- Durable & Patient Ergonomic Design
- Lightweight & Outperforming Quality Engineering
- Innovative & Convenient Mask Setup Positioning

**MR  
PROTON**

**• Included components**

**PROMISE™ Baseplate (MR)**  
R680-AGF



**PROMISE™ Head Mask**  
RG-PH-2436C (Green)

**PROMISE™ Couch Insert for  
Treatment Delivery**  
R680-ACFP



**PROMISE™ Head AccuCushion**  
RG-550-CEL2 (Green)





MR Patient  
Positioning Devices

Klarity ATHENA



Klarity ATHENA™ Supine Breast System  
R610-PDKF



MR  
ATHENA

### Klarity ATHENA™ Supine Breast System

- Well accommodated in 70cm bore
- More reliable with max load limit 200 kgs
- Enables the easy use with vacuum bags, head mask or head & thorax mask for better personalized immobilization



Klarity ATHENA™ Breast Wing Board (MR Version)  
R632-NTKF





**MR Conditional up to 3.0T**

### Klarity MR Vacuum Bags

- The slip-proof surface increases the bag stability when fixing on the overlay
- The bags can use with many kinds of immobilization devices
- Customization is available

**MR  
Vacuum Bags**

## Ordering Information

Model	Description	Specification	Compatibility
MR7202-9N	Vaccum Bag for Head & Shoulder Support	65.2X68.5 cm	Universal
MR7206-16N	Vaccum Bag for Head & Shoulder Support	98.5X70 cm	Universal
MR7244-21NR-O	Vaccum Bag for Head & Shoulder Support	135.2X60 cm	For R612 & R605 Baseplate
MR7249-21NR-O	Vaccum Bag for Head & Shoulder Support	134.7X70 cm	For R612 & R605 Baseplate
MR7504-28N	Vaccum Bag for Body Support	100X70 cm	Universal
MR7507-36N	Vaccum Bag for Body Support	120X80 cm	Universal
MR7580-36N	Vaccum Bag for Body Support	140X70 cm	Universal
MR75103-77N	Vaccum Bag for Body Support	200X80 cm	Universal
MR7644-40N	Vaccum Bag for Body Support	142X92 cm	Universal
MR7641-25NB4	Vaccum Bag for Body Support	120X80 cm	For SBRT baseplate only
MR7704-5NT-O	Vaccum Bag for Limbs Support	80X28 cm	For SBRT baseplate only



## MR Baseplates



**MR Optek™ System**  
R630-3SGF



**MR Optek™ Demi-Baseplate**  
R630-SHGF  
Compatible with Elekta Elekta  
HexaPOD™ System




MR Conditional up to 3.0T

**Prone Pelvis System**  
R612-AE4




**Klarity Leg Positioner**  
R634-LGF

## MR Accessories



**Foam Cushion**  
R516-178L  
139x178x381mm



### MR Indexing Bars

**RE-17MGF**

Glass Fiber Indexing Bar for Elekta Unity MR Linac

**RE-630ESGF**

Glass Fiber Indexing Bar for Elekta iBeam, TomoTherapy & United Imaging uRT-linac 506c

**RE-7BGF**

Glass Fiber Indexing Bar for Varian Truebeam

**RE-7SGF**

Glass Fiber Indexing Bar for Varian Truebeam  
Glass Fiber Indexing Bar for Elekta Unity MR Linac



**MR Head Support Set of 6**  
R507-GF



KLARTY CHINA  
Klarity Medical & Equipment (GZ) Co., Ltd.  
No.14, 3rd Street Shawan  
GETDD Guangzhou 510730,China  
Tel: +86-20-8222 6380  
Fax: +86-20-8208 8552  
Sales@klarity-medical.com

Version Number: KL-CY-20220715-A0



LinkedIn



YouTube