

Klarity[®]



Asia & Worldwide
Klarity Medical & Equipment Co., Ltd
No. 14, 3rd Street Shawan
GETDD, 510730, Guangzhou
Guangdong, China
Tel +86-20-8222 6380
Fax +86-20-8208 8552
sales@klarity-medical.com
www.klarity-medical.com/en

North America
Klarity Medical Products USA
600 Industrial Parkway
Heath, Ohio 43056
Tel (877) 599-8905 (toll free)
Tel (740) 788-8107
Fax (740) 788-8109
hello@klaritymedical.com
www.klaritymedical.com




Klarity on LinkedIn

Klarity UNO[™] All-In-One System



All-In-One System



Klarity UNO™ All-In-One Carbon Fibre Baseplate

Compatible with most positioning devices



Klarity UNO™ All-In-One Carbon Fibre Baseplate Code: R616-SCFA

MR Klarity UNO™ All-In-One MR Compatible Baseplate Code: R616-SKFA

Compatible with small bore CT, MR, and Elekta Unity



Klarity UNO™ All-In-One Carbon Fibre Baseplate Code: R616-NCF

MR Klarity UNO™ All-In-One MR Compatible Baseplate Code: R616-NKF

Compatible with Elekta HexaPod




Klarity UNO™ All-In-One Carbon Fibre Baseplate Code: R616-SCF1A

MR Klarity UNO™ All-In-One MR Compatible Baseplate Code: R616-SKF1A

All-In-One Baseplates with Standard Accessories

- > Improved Efficiency: Streamlined workflow through full body positioning
- > Lightweight & Strong: Supports up to 200kg



Klarity UNO™ Head & Neck Solution



> Compatible with S-type Head-only and Head & Shoulder Masks, as well as Klarity Cushions™ and Vacuum Bags




Head & Neck Solution Accessories



- 



Head Masks


- 

Head & Shoulder Masks


- 

Klarity Cushions™

- 

Head & Shoulder Klarity Cushions™

- 

Head & Shoulder Headrest R507-T

- 

Cushion Molding Tool RSX-K612SA



Klarity UNO™ All-In-One Carbon Fibre Baseplate

Klarity UNO™

Head & Thorax and Hip & Pelvis Solutions with Overhead Arm Positioner

- > Multiple positions for patient immobilization
- > Excellent comfort and precise reproducibility



Head & Thorax and Hip & Pelvis Solution Accessories



Overhead Arm Positioner
R616-WCFA



Leg Positioner Cushion
R616-L1



Foot Holder
R616-L3



Head Supports



Head & Thorax Masks



Thorax Masks




Pelvis Masks



Klarity Cushions™



Cushion Indexer
RSX-TZ-616A



Klarity UNO™ All-In-One
Carbon Fibre Baseplate

Klarity UNO™





Head & Breast Solution

> Offers flexible solutions for better customized positioning



Klarity UNO™ All-In-One Carbon Fibre Baseplate

Head & Breast Solution Accessories



Klarity UNO™ Head & Torso Solution

> Provides comprehensive immobilization for the whole brain and spinal cord



Klarity UNO™ All-In-One
Carbon Fibre Baseplate

Craniospinal Irradiation Whole Body Fixation Accessories



Head & Shoulder Masks



Head Supports



Head Supports



Long Torso Masks

Klarity UNO™

SRS Solution

- > Frameless, non-invasive
- > Sub-millimeter accuracy
- > Suitable for 4Pi SRS immobilization



SRS Solution Accessories



SRS Adaptor
R616-ACF

MR




SRS Mask Set, Basic
RBSRS-F01

MR



SRS Mask Set, Premium
RBSRS-F02

MR



Klarity UNO™ All-In-One
Carbon Fibre Baseplate

Klarity UNO™ SBRT Solution

> User-friendly with multiple abdominal compression options



Klarity UNO™ All-In-One
Carbon Fibre Baseplate

SBRT Solution Accessories



SBRT Solution Accessories
R616-BT1



Compression Belt
R616-AIR



Thorax Mask



SBRT Abdominal Compression Solution

- > Easy-to-use Belly Bridge with adjustable multi-directional support
- > Enhances treatment precision by measuring and reducing motion
- > Specially designed curved shape for seamless integration with O-ring radiation therapy equipment



Flexible Compression Solution
Better patient comfort



Abdominal Compression Mask Solution
Improves intra-fractional stability during treatment



Compression Belt Solution
Reduces attenuation through abdominal compression




Klarity UNO™

Prone Pelvis Solution

- > Designed for better patient comfort
- > Enhanced accuracy with additional positioning accessories

The Prone Breast Solution can be used with the Klarity UNO baseplate or it can be used on its own.

Prone Pelvis Solution Accessories



Klarity UNO™ All-In-One Carbon Fibre Baseplate

Klarity UNO™


Prone Breast Solution

- > Designed for better patient comfort
- > Flexible positioning for right or left breast



The Prone Breast Solution can be used with the Klarity UNO baseplate or it can be used on its own.

Prone Breast Solution Accessories



Prone Breast Cushion Set
R616-BE5




Thorax Masks



Breast Klarity Cushions™



Prone Breast Vacuum Bag
R7652-28NL-O



Klarity UNO™ All-In-One
Carbon Fibre Baseplate

R616 Quick Release System



Our quick release system is designed to meet the needs of patient comfort during treatment. Doctors can effortlessly unlock the mask with a single press of the button, empowering patients to move freely and comfortably.

R616 Quick Release System Accessories




Klarity UNO™ All-In-One Carbon Fibre Baseplate



R616 Quick Release System Accessories R616-SCFR1-A



Head Masks



Head & Shoulder Masks




Klarity Cushions™



Head & Shoulder Klarity Cushions™

Uno All-In-One Storage Cabinet Solution




Clinical Adoption



Consumables and Accessories



S-type Head Masks
Compatible with all standard S-Type Head Masks



S-type Head & Shoulder Masks
Compatible with all standard S-Type Head Masks



Head & Thorax Masks
RLG322-S42436 (LiteGreen™)
R322-S42436 (White)
with 36% perforation, 2.4mm



Head & Thorax Masks
RLG322-S424C-L (LiteGreen™)
R322-S424C-L (White)
Cutout on Left side




Head & Thorax Masks
RLG322-S424C-R (LiteGreen™)
R322-S424C-R (White)
Cutout on Right side



Thorax Masks
RLG902-2436TL (LiteGreen™)
R902-2436TL (White)
with 36% perforation, 2.4mm



Pelvis Masks
RLG902-2436BA (LiteGreen™)
R902-2436BA (White)
with 36% perforation, 2.4mm



Long Torso Masks
RLG902-2436LA (LiteGreen™)
R902-2436LA (White)
with 36% perforation, 2.4mm



Head & Shoulder Vacuum Bag
R7202-9NLR
625 X 582 mm, 9 Liters fill



Body Vacuum Bag
R7569-12NLR-O
945 X 530 mm, 12 Liters fill



Full Body Vacuum Bag
R7231-15NLR-O
1347 X 520 mm, 15 Liters fill
Full Body Vacuum Bag (SBRT)
R7231-15NLR-O
1350 X 700 mm, 21 Liters fill



Prone Breast Vacuum Bag
R7652-28NL-O
930 X 810 mm, 28 Liters fill



Prone Pelvis Vacuum Bag
R7653-20NL-O
900 X 850 mm, 20 Liters fill



Vacuum Bag for Body support
R7641-25NLB5
1200 X 800 mm, 25 Liters fill



Vacuum Bag for Knee Support
R7136-4NLT-O
800 X 280 mm, 4 Liters fill



Prone Breast Vacuum Bag
MR7653-20NL-O
900 X 850 mm, 20 Liters fill



Body Vacuum Bags
MR7677-25N
25 Liters fill



Klarity Cushions™
Compatible with all Klarity Cushions™ and AccuCushions®



Head & Shoulder Klarity Cushions™
Compatible with all H&S Klarity Cushions™ and AccuCushions®



Pelvis Klarity Cushions™
Compatible with all H&S Klarity Cushions™ and AccuCushions®



Breast Klarity Cushions™
Compatible with all H&S Klarity Cushions™ and AccuCushions®



Body Klarity Cushions™
Compatible with all H&S Klarity Cushions™ and AccuCushions®



Head & Shoulder Headrest
R507-T (276 X 424 X 56 mm)



Cushion Molding Tool
RSX-K612SA (450 X 478 X 60 mm)



Mask Adaptor
R616-NTCF-3



Cushion Indexer
RSX-TZ-616A



R616 Respiratory Belt
R616-AIR